3mach®-Kappenadapter

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3mach® cap adapters

Use: Accessories for 3mach® flushing modules in the area of instrument preparation.

3mach® cap adapters (Luer lock, plug-in and tube adapters, with and without a stainless steel thread) are designed for cleaning medical instruments during the cleaning process of a cleaning and disinfection machine. Ensure that the instrument diameters or plug-in adapters are correctly inserted in the intended silicone cap apertures. 3mach® silicone caps with valve function can seal themselves and have an internal flap for this purpose. In addition, manufacturer specifications on instrument preparation must be obtained from the instrument manufacturer.

Caution: All 3mach® cap adapters must be cleaned prior to initial use. Screw-fitted cap adapters may loosen due to high temperature fluctuations. Ensure that connection does not cause the Luer lock adapters to tilt relative to each other to enable their release by hand. Silicone caps are damaged by sharp-edged instruments and should therefore not be introduced. Never use 3mach® cap adapters on people.

Note: 3mach® cap adapters require 3mach® silicone caps to clean medical lumen instruments. Residual water in adapters and tubes must be removed. Brittle or defective silicone caps must be replaced.

Preparation: Suitable for alkaline cleaning processes in the area of mechanical cleaning, thermal disinfection and drying (up to 120 °C) and use in medical technology during fractionated vacuum processes (134 °C). For further instructions for use or types of use, please contact us by email at info@3mach.de.

Liability: 3mach GmbH bears no liability for injury or damage resulting from the use of 3mach® products. This also applies to individually converted and modified 3mach® products and their misuse.

Manufacturer: 3mach GmbH, Daimlerstrasse 25, 72414 Rangendingen, Germany
Tel.: 0049 (0) 74 71-98 41 18, Fax: 0049 (0) 74 71-98 41 19, www.3mach.de

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Hersteller: Firma 3mach GmbH, Daimlerstraße 25, 72414 Rangendingen, Germany